

Viscosity

OFITE DIGITAL STOPWATCH #155-25

To properly perform the Marsh Funnel viscosity test, the OFITE electronic digital electronic stopwatch may be used. The quartz crystal timer counts up to 30 minutes in .01 seconds and counts in seconds if over 30 minutes. Each stopwatch is calibrated against the National Institute of Standards and a test certificate is supplied with each timer. The calibration complies fully with the requirements of ISO 9000 Certification.

Size: 2.4" x 2" x 0.6" (6 x 5 x 2 cm)
Weight: 1.5 oz (42.5 grams)



*Did you
know?*

When you send something to us for repair, we know you are without a critical piece of equipment. That's why we make your repairs first priority (even equipment from other manufacturers!).

Also Available:

- #155-26 Basic Stopwatch, Bodytronics
- #110-30 Measuring Cup, Stainless Steel, 500 mL
- #110-40 Measuring Cup, Stainless Steel, 1,000 mL
- #110-50 Measuring Cup, Stainless Steel, 2,000 mL



Stainless Steel Measuring Cups