

# HTHP Filtration

## HTHP SINGLE-END CELL ASSEMBLY FOR CEMENT TESTING, 175 ML, 2,000 PSI #170-45-1

**Size:** 3.25" x 3.25" x 4.5" (8 x 8 x 11 cm)  
**Weight:** 7 lb 14 oz (3.6 kg)

### Components:

- #170-13 O-ring for Cell, Buna N
- #170-16 Valve Stem
- #170-17 O-ring for Valve Stem, Viton®
- #170-18 Screen, 325-Mesh with 60-Mesh Backup, Detachable
- #170-24 End Cap for Detachable Screen, 2,000 PSI
- #170-26 Cap Locking Screw, Stainless Steel
- #170-27 Allen Wrench for Cap Locking Screw, 5/32"
- #170-45-2 Cell Body, Single-End, 2,000 PSI, 316 Stainless Steel



#170-45-1 Cell Assembly for Cement Testing

## HTHP SINGLE-END CELL ASSEMBLY FOR CERAMIC DISKS, 175 ML, 2,000 PSI #170-48

**Size:** 3.25" x 3.25" x 4.5" (8 x 8 x 11 cm)  
**Weight:** 7 lb 14 oz (3.6 kg)

### Components:

- #170-13 O-ring for Cell, Buna N
- #170-16 Valve Stem
- #170-17 O-ring for Valve Stem, Viton®
- #170-26-1 Hardened Cap Locking Screw
- #170-27 Allen Wrench for Cap Locking Screw, 5/32"
- #170-49 Cell Body, 175 mL
- #170-69 Scribed Cell Cap for Ceramic Disk, 2,000 PSI



#170-48 Cell Assembly for Ceramic Disks, 175 mL

## HTHP DOUBLE-END CELL ASSEMBLY FOR CERAMIC DISKS, 175 ML, 303 STAINLESS STEEL, 2,000 PSI #170-46

**Size:** 3.25" x 3.25" x 4.5" (8 x 8 x 11 cm)  
**Weight:** 7 lb 14 oz (3.6 kg)

### Components:

- #170-13 O-ring for Cell, Buna N
- #170-16 Valve Stem
- #170-17 O-ring for Valve Stem, Viton®
- #170-26-1 Hardened Cap Locking Screw
- #170-27 Allen Wrench for Cap Locking Screw, 5/32"
- #170-47 Cell Body for Ceramic Disk, 175 mL
- #170-69 Cell Cap for Ceramic Disk, Scribed, 2,000 PSI
- #170-72 Spacer for Filter Paper, 1/4", Stainless Steel
- #170-77 O-ring for Stainless Steel Spacer
- #171-21 Cell Cap with Screen, 2,000 PSI



#170-46 Cell Assembly for Ceramic Disks, 175 mL